




Meditech

Meditech Ortho India Company



Locking Plate System

Manufacturer & Exporter
Orthopaedic Implants & Instruments

FDA
CERTIFIED

CE
CERTIFIED

ISO 9001:2008
ISO 13484:2012

CONTENT

- **Mini Fragment** ----- **01 - 26**
- **Small Fragment** ----- **27 - 55**
- **Large Fragment** ----- **56 - 77**
- **Wise Locked Screw** ----- **78 - 86**

2.7mm Distal Radius Volar Plate 8Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM
(Prod. Code)


SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 201.02	MS 201.02	2 Hole (L/R)
MT 201.03	MS 201.03	3 Hole (L/R)
MT 201.04	MS 201.04	4 Hole (L/R)
MT 201.05	MS 201.05	5 Hole (L/R)
MT 201.06	MS 201.06	6 Hole (L/R)
MT 201.07	MS 201.07	7 Hole (L/R)



Meditech



2.7mm Distal Radius Volar Plate 9Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM
(Prod. Code)


SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 202.02	MS 202.02	2 Hole (L/R)
MT 202.03	MS 202.03	3 Hole (L/R)
MT 202.04	MS 202.04	4 Hole (L/R)
MT 202.05	MS 202.05	5 Hole (L/R)
MT 202.06	MS 202.06	6 Hole (L/R)
MT 202.07	MS 202.07	7 Hole (L/R)



Meditech



2.7mm Distal Radius Volar Plate Extra Articular 4Hole Head

Primary Use : Distal Radius Fracture


Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM
(Prod. Code)

SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 203.02	MS 203.02	2 Hole (L/R)
MT 203.03	MS 203.03	3 Hole (L/R)
MT 203.04	MS 203.04	4 Hole (L/R)
MT 203.05	MS 203.05	5 Hole (L/R)
MT 203.06	MS 203.06	6 Hole (L/R)
MT 203.07	MS 203.07	7 Hole (L/R)




2.7mm Distal Radius Volar Plate Extra Articular Long 4Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 204.05	MS 204.05	05 Hole (L/R)
MT 204.06	MS 204.06	06 Hole (L/R)
MT 204.07	MS 204.07	07 Hole (L/R)
MT 204.08	MS 204.08	08 Hole (L/R)
MT 204.09	MS 204.09	09 Hole (L/R)
MT 204.10	MS 204.10	10 Hole (L/R)
MT 204.11	MS 204.11	11 Hole (L/R)
MT 204.12	MS 204.12	12 Hole (L/R)



Meditech



2.7mm Distal Radius Volar Plate Extra Articular 5Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 205.02	MS 205.02	2 Hole (L/R)
MT 205.03	MS 205.03	3 Hole (L/R)
MT 205.04	MS 205.04	4 Hole (L/R)
MT 205.05	MS 205.05	5 Hole (L/R)
MT 205.06	MS 205.06	6 Hole (L/R)
MT 205.07	MS 205.07	7 Hole (L/R)



Meditech



2.7mm Distal Radius Volar Buttress Plate 5Hole Head

Primary Use : Distal Radius Fracture


Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM
(Prod. Code)

SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 206.02	MS 206.02	2 Hole (L/R)
MT 206.03	MS 206.03	3 Hole (L/R)
MT 206.04	MS 206.04	4 Hole (L/R)
MT 206.05	MS 206.05	5 Hole (L/R)
MT 206.06	MS 206.06	6 Hole (L/R)
MT 206.07	MS 206.07	7 Hole (L/R)



Meditech



2.7mm Distal Radius Angular Plate 5Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM
(Prod. Code)

SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 207.02

MS 207.02

2 Hole (L/R)

MT 207.03

MS 207.03

3 Hole (L/R)

MT 207.04

MS 207.04

4 Hole (L/R)

MT 207.05

MS 207.05

5 Hole (L/R)

MT 207.06

MS 207.06

6 Hole (L/R)

MT 207.07

MS 207.07

7 Hole (L/R)



Meditech




2.7mm Distal Radius T-Wrist Plate

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 208.03	MS 208.03	03 Hole (L/R)
MT 208.04	MS 208.04	04 Hole (L/R)
MT 208.05	MS 208.05	05 Hole (L/R)
MT 208.06	MS 208.06	06 Hole (L/R)
MT 208.07	MS 208.07	07 Hole (L/R)
MT 208.08	MS 208.08	08 Hole (L/R)
MT 208.09	MS 208.09	09 Hole (L/R)
MT 208.10	MS 208.10	10 Hole (L/R)
MT 208.11	MS 208.11	11 Hole (L/R)
MT 208.12	MS 208.12	12 Hole (L/R)



Meditech




2.7mm Distal Radius Volar Plate Variable Angle


Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 209.02	MS 209.02	02 Hole (L/R)
MT 209.03	MS 209.03	03 Hole (L/R)
MT 209.04	MS 209.04	04 Hole (L/R)
MT 209.05	MS 209.05	05 Hole (L/R)
MT 209.06	MS 209.06	06 Hole (L/R)
MT 209.07	MS 209.07	07 Hole (L/R)
MT 209.08	MS 209.08	08 Hole (L/R)



Meditech



2.7mm Distal Radius L-Dorsal Locking Plate 3Hole Head


Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 210.02	MS 210.02	02 Hole (L/R)
MT 210.03	MS 210.03	03 Hole (L/R)
MT 210.04	MS 210.04	04 Hole (L/R)
MT 210.05	MS 210.05	05 Hole (L/R)
MT 210.06	MS 210.06	06 Hole (L/R)
MT 210.07	MS 210.07	07 Hole (L/R)
MT 210.08	MS 210.08	08 Hole (L/R)



Meditech



2.7mm Distal Radius Oblique L-Dorsal Locking Plate 3Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM
(Prod. Code)


SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 211.02	MS 211.02	02 Hole (L/R)
MT 211.03	MS 211.03	03 Hole (L/R)
MT 211.04	MS 211.04	04 Hole (L/R)
MT 211.05	MS 211.05	05 Hole (L/R)
MT 211.06	MS 211.06	06 Hole (L/R)
MT 211.07	MS 211.07	07 Hole (L/R)
MT 211.08	MS 211.08	08 Hole (L/R)



Meditech



2.7mm Distal Radius L-Dorsal Locking Plate 2Hole Head

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 212.02	MS 212.02	02 Hole (L/R)
MT 212.03	MS 212.03	03 Hole (L/R)
MT 212.04	MS 212.04	04 Hole (L/R)
MT 212.05	MS 212.05	05 Hole (L/R)
MT 212.06	MS 212.06	06 Hole (L/R)
MT 212.07	MS 212.07	07 Hole (L/R)
MT 212.08	MS 212.08	08 Hole (L/R)



Meditech



2.7mm Distal Radius Oblique L-Dorsal Locking Plate 2Hole Head


Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 213.02	MS 213.02	02 Hole (L/R)
MT 213.03	MS 213.03	03 Hole (L/R)
MT 213.04	MS 213.04	04 Hole (L/R)
MT 213.05	MS 213.05	05 Hole (L/R)
MT 213.06	MS 213.06	06 Hole (L/R)
MT 213.07	MS 213.07	07 Hole (L/R)
MT 213.08	MS 213.08	08 Hole (L/R)



Meditech




2.7mm Distal Radius Dorsal Locking Plate Straight


Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 214.02	MS 214.02	02 Hole
MT 214.03	MS 214.03	03 Hole
MT 214.04	MS 214.04	04 Hole
MT 214.05	MS 214.05	05 Hole
MT 214.06	MS 214.06	06 Hole
MT 214.07	MS 214.07	07 Hole
MT 214.08	MS 214.08	08 Hole



Meditech



2.7mm Small T-Locking Plate

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 215.02	MS 215.02	02 Hole
MT 215.03	MS 215.03	03 Hole
MT 215.04	MS 215.04	04 Hole
MT 215.05	MS 215.05	05 Hole
MT 215.06	MS 215.06	06 Hole
MT 215.07	MS 215.07	07 Hole
MT 215.08	MS 215.08	08 Hole



Meditech




2.7mm Small T-Adaption Locking Plate

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 216.02	MS 216.02	02 Hole
MT 216.03	MS 216.03	03 Hole
MT 216.04	MS 216.04	04 Hole
MT 216.05	MS 216.05	05 Hole
MT 216.06	MS 216.06	06 Hole
MT 216.07	MS 216.07	07 Hole
MT 216.08	MS 216.08	08 Hole



Meditech



2.7mm Small T-Adaption Locking Plate Round Hole

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

17

TITANIUM
(Prod. Code)**SS. 316L**
Prod. Code**SPECIFICATION**
(Holes)

MT 217.02	MS 217.02	02 Hole
MT 217.03	MS 217.03	03 Hole
MT 217.04	MS 217.04	04 Hole
MT 217.05	MS 217.05	05 Hole
MT 217.06	MS 217.06	06 Hole
MT 217.07	MS 217.07	07 Hole
MT 217.08	MS 217.08	08 Hole

Meditech



2.0mm Small Struck Locking Plate

Primary Use : Distal Radius Fracture

Matching Screw : 2.0mm Locking and 2.0mm Cortical Screw

18

TITANIUM
(Prod. Code)**SS. 316L**
Prod. Code**SPECIFICATION**
(Holes)

MT 218.02	MS 218.02	02 Hole
MT 218.03	MS 218.03	03 Hole
MT 218.04	MS 218.04	04 Hole
MT 218.05	MS 218.05	05 Hole
MT 218.06	MS 218.06	06 Hole
MT 218.07	MS 218.07	07 Hole
MT 218.08	MS 218.08	08 Hole

Meditech



2.7mm Mini Condylar Locking Plate Round Hole


Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 219.02	MS 219.02	02 Hole (L/R)
MT 219.03	MS 219.03	03 Hole (L/R)
MT 219.04	MS 219.04	04 Hole (L/R)
MT 219.05	MS 219.05	05 Hole (L/R)
MT 219.06	MS 219.06	06 Hole (L/R)
MT 219.07	MS 219.07	07 Hole (L/R)
MT 219.08	MS 219.08	08 Hole (L/R)



Meditech



2.7mm Mini Locking Plate

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

20

TITANIUM
(Prod. Code)**SS. 316L**
Prod. Code**SPECIFICATION**
(Holes)

MT 220.02

MS 220.02

02 Hole

MT 220.03

MS 220.03

03 Hole

MT 220.04

MS 220.04

04 Hole

MT 220.05

MS 220.05

05 Hole

MT 220.06

MS 220.06

06 Hole

MT 220.07

MS 220.07


07 Hole

MT 220.08

MS 220.08

08 Hole

Meditech



2.7mm Adaption Locking Plate

Primary Use : Distal Radius Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 221.02	MS 221.02	02 Hole
MT 221.03	MS 221.03	03 Hole
MT 221.04	MS 221.04	04 Hole
MT 221.05	MS 221.05	05 Hole
MT 221.06	MS 221.06	06 Hole
MT 221.07	MS 221.07	07 Hole
MT 221.08	MS 221.08	08 Hole



Meditech




2.7mm Proximal Radial Head Locking Plate

Primary Use : Radial Head Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 222.02	MS 222.02	02 Hole (L/R)
MT 222.03	MS 222.03	03 Hole (L/R)
MT 222.04	MS 222.04	04 Hole (L/R)
MT 222.05	MS 222.05	05 Hole (L/R)
MT 222.06	MS 222.06	06 Hole (L/R)
MT 222.07	MS 222.07	07 Hole (L/R)
MT 222.08	MS 222.08	08 Hole (L/R)




2.7mm Proximal Radial Head Buttress Locking Plate

Primary Use : Radial Head Fracture

Matching Screw : 2.7mm Locking and 2.7mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 223.02	MS 223.02	02 Hole (L/R)
MT 223.03	MS 223.03	03 Hole (L/R)
MT 223.04	MS 223.04	04 Hole (L/R)
MT 223.05	MS 223.05	05 Hole (L/R)
MT 223.06	MS 223.06	06 Hole (L/R)
MT 223.07	MS 223.07	07 Hole (L/R)
MT 223.08	MS 223.08	08 Hole (L/R)



Meditech



2.0mm Mini H-Locking Plate

Primary Use : Small Osteotomies Fracture

Matching Screw : 2.0mm Locking Screw

24

TITANIUM
(Prod. Code)**SS. 316L**
Prod. Code**SPECIFICATION**
(Holes)

MT 224.04

MS 224.04

04 Hole



Meditech




3.5mm Small T-Locking Plate


Primary Use : Distal Radius Fracture

Matching Screw : 3.5mm Locking and 3.5mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 225.02	MS 225.02	02 Hole
MT 225.03	MS 225.03	03 Hole
MT 225.04	MS 225.04	04 Hole
MT 225.05	MS 225.05	05 Hole
MT 225.06	MS 225.06	06 Hole
MT 225.07	MS 225.07	07 Hole
MT 225.08	MS 225.08	08 Hole



Meditech



3.5mm Small T-Oblique Locking Plate

Primary Use : Distal Radius Fracture


Matching Screw : 3.5mm Locking and 3.5mm Cortical Screw

TITANIUM
(Prod. Code)


SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 226.02	MS 226.02	02 Hole (L/R)
MT 226.03	MS 226.03	03 Hole (L/R)
MT 226.04	MS 226.04	04 Hole (L/R)
MT 226.05	MS 226.05	05 Hole (L/R)
MT 226.06	MS 226.06	06 Hole (L/R)
MT 226.07	MS 226.07	07 Hole (L/R)
MT 226.08	MS 226.08	08 Hole (L/R)



Meditech




3.5mm Small Locking Plate

Primary Use : Both Bone Fracture

Matching Screw : 3.5mm Locking and 3.5mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 227.04	MS 227.04	04 Hole
MT 227.05	MS 227.05	05 Hole
MT 227.06	MS 227.06	06 Hole
MT 227.07	MS 227.07	07 Hole
MT 227.08	MS 227.08	08 Hole
MT 227.09	MS 227.09	09 Hole
MT 227.10	MS 227.10	10 Hole
MT 227.11	MS 227.11	11 Hole
MT 227.12	MS 227.12	12 Hole

Meditech




1/3 Tubular Locking Plate

Primary Use : Both Bone and Distal Fibular Fracture

Matching Screw : 3.5mm Locking and 3.5mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 228.04	MS 228.04	04 Hole
MT 228.05	MS 228.05	05 Hole
MT 228.06	MS 228.06	06 Hole
MT 228.07	MS 228.07	07 Hole
MT 228.08	MS 228.08	08 Hole
MT 228.09	MS 228.09	09 Hole
MT 228.10	MS 228.10	10 Hole
MT 228.11	MS 228.11	11 Hole
MT 228.12	MS 228.12	12 Hole

Meditech




3.5mm Recon Locking Plate Low Profile

Primary Use : Pelvic and Clavicle Fracture

Matching Screw : 3.5mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 229.04	MS 229.04	04 Hole
MT 229.05	MS 229.05	05 Hole
MT 229.06	MS 229.06	06 Hole
MT 229.07	MS 229.07	07 Hole
MT 229.08	MS 229.08	08 Hole
MT 229.09	MS 229.09	09 Hole
MT 229.10	MS 229.10	10 Hole
MT 229.11	MS 229.11	11 Hole
MT 229.12	MS 229.12	12 Hole
MT 229.13	MS 229.13	13 Hole
MT 229.14	MS 229.14	14 Hole
MT 229.15	MS 229.15	15 Hole



Meditech




3.5mm Recon Locking Plate

Primary Use : Pelvic and Clavicle Fracture

Matching Screw : 3.5mm Locking Screw and 3.5mm Cortical Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 230.04	MS 230.04	04 Hole
MT 230.05	MS 230.05	05 Hole
MT 230.06	MS 230.06	06 Hole
MT 230.07	MS 230.07	07 Hole
MT 230.08	MS 230.08	08 Hole
MT 230.09	MS 230.09	09 Hole
MT 230.10	MS 230.10	10 Hole
MT 230.11	MS 230.11	11 Hole
MT 230.12	MS 230.12	12 Hole
MT 230.13	MS 230.13	13 Hole
MT 230.14	MS 230.14	14 Hole
MT 230.15	MS 230.15	15 Hole




3.5mm Recon Curved Locking Plate

Primary Use : Pelvic Fracture

Matching Screw : 3.5mm Locking Screw and 3.5mm Cortical Screw

31

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 231.04	MS 231.04	04 Hole
MT 231.05	MS 231.05	05 Hole
MT 231.06	MS 231.06	06 Hole
MT 231.07	MS 231.07	07 Hole
MT 231.08	MS 231.08	08 Hole
MT 231.09	MS 231.09	09 Hole
MT 231.10	MS 231.10	10 Hole
MT 231.11	MS 231.11	11 Hole
MT 231.12	MS 231.12	12 Hole
MT 231.13	MS 231.13	13 Hole
MT 231.14	MS 231.14	14 Hole
MT 231.15	MS 231.15	15 Hole
MT 231.16	MS 231.16	16 Hole




Olecranon Locking Plate


Primary Use : Olecranon Fracture

Matching Screw : 3.5mm Locking, 3.5mm Cortical and 4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 232.03	MS 232.03	03 Hole (L/R)
MT 232.04	MS 232.04	04 Hole (L/R)
MT 232.05	MS 232.05	05 Hole (L/R)
MT 232.06	MS 232.06	06 Hole (L/R)
MT 232.07	MS 232.07	07 Hole (L/R)
MT 232.08	MS 232.08	08 Hole (L/R)
MT 232.09	MS 232.09	09 Hole (L/R)
MT 232.10	MS 232.10	10 Hole (L/R)
MT 232.11	MS 232.11	11 Hole (L/R)
MT 232.12	MS 232.12	12 Hole (L/R)



Meditech



Distal Medial Humerus Locking Plate

Primary Use : Distal humerus Fracture

Matching Screw : 3.5mm Locking, 3.5mm Cortical and 2.7mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 233.03	MS 233.03	03 Hole (L/R)
MT 233.04	MS 233.04	04 Hole (L/R)
MT 233.05	MS 233.05	05 Hole (L/R)
MT 233.06	MS 233.06	06 Hole (L/R)
MT 233.07	MS 233.07	07 Hole (L/R)
MT 233.08	MS 233.08	08 Hole (L/R)
MT 233.09	MS 233.09	09 Hole (L/R)
MT 233.10	MS 233.10	10 Hole (L/R)
MT 233.11	MS 233.11	11 Hole (L/R)
MT 233.12	MS 233.12	12 Hole (L/R)

Meditech




Dorsolateral Distal Humerus Locking Plate with Support


Primary Use : Distal humerus Fracture

Matching Screw : 3.5mm Locking, 3.5mm Cortical and 2.7mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 234.03	MS 234.03	03 Hole (L/R)
MT 234.04	MS 234.04	04 Hole (L/R)
MT 234.05	MS 234.05	05 Hole (L/R)
MT 234.06	MS 234.06	06 Hole (L/R)
MT 234.07	MS 234.07	07 Hole (L/R)
MT 234.08	MS 234.08	08 Hole (L/R)
MT 234.09	MS 234.09	09 Hole (L/R)
MT 234.10	MS 234.10	10 Hole (L/R)
MT 234.11	MS 234.11	11 Hole (L/R)
MT 234.12	MS 234.12	12 Hole (L/R)



Meditech




Dorsolateral Distal Humerus Locking Plate without Support


Primary Use : Distal humerus Fracture

Matching Screw : 3.5mm Locking, 3.5mm Cortical and 2.7mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 235.03	MS 235.03	03 Hole (L/R)
MT 235.04	MS 235.04	04 Hole (L/R)
MT 235.05	MS 235.05	05 Hole (L/R)
MT 235.06	MS 235.06	06 Hole (L/R)
MT 235.07	MS 235.07	07 Hole (L/R)
MT 235.08	MS 235.08	08 Hole (L/R)
MT 235.09	MS 235.09	09 Hole (L/R)
MT 235.10	MS 235.10	10 Hole (L/R)
MT 235.11	MS 235.11	11 Hole (L/R)
MT 235.12	MS 235.12	12 Hole (L/R)



Meditech




Extra Articular Distal Humerus Locking Plate

Primary Use : Extra Articular Distal humerus Fracture

Matching Screw : 3.5mm Locking, 3.5mm Cortical and 4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 236.03	MS 236.03	03 Hole (L/R)
MT 236.04	MS 236.04	04 Hole (L/R)
MT 236.05	MS 236.05	05 Hole (L/R)
MT 236.06	MS 236.06	06 Hole (L/R)
MT 236.07	MS 236.07	07 Hole (L/R)
MT 236.08	MS 236.08	08 Hole (L/R)
MT 236.09	MS 236.09	09 Hole (L/R)
MT 236.10	MS 236.10	10 Hole (L/R)
MT 236.11	MS 236.11	11 Hole (L/R)
MT 236.12	MS 236.12	12 Hole (L/R)



Meditech




Metaphysel Distal Humerus Locking Plate

Primary Use : Condylar Distal humerus Fracture

Matching Screw : 3.5mm Locking, 3.5mm Cortical and 4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 237.03	MS 237.03	03 Hole (L/R)
MT 237.04	MS 237.04	04 Hole (L/R)
MT 237.05	MS 237.05	05 Hole (L/R)
MT 237.06	MS 237.06	06 Hole (L/R)
MT 237.07	MS 237.07	07 Hole (L/R)
MT 237.08	MS 237.08	08 Hole (L/R)
MT 237.09	MS 237.09	09 Hole (L/R)
MT 237.10	MS 237.10	10 Hole (L/R)
MT 237.11	MS 237.11	11 Hole (L/R)
MT 237.12	MS 237.12	12 Hole (L/R)



Meditech




Metaphysel Distal Tibia Locking Plate**Primary Use :** Lower and Tibia Fracture**Matching Screw :** 3.5mm Cortical and 4.0mm Cancellous Locking Screw

4.5mm Cortical and 5.0mm Locking Screw

38

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 238.03	MS 238.03	03 Hole (L/R)
MT 238.04	MS 238.04	04 Hole (L/R)
MT 238.05	MS 238.05	05 Hole (L/R)
MT 238.06	MS 238.06	06 Hole (L/R)
MT 238.07	MS 238.07	07 Hole (L/R)
MT 238.08	MS 238.08	08 Hole (L/R)
MT 238.09	MS 238.09	09 Hole (L/R)
MT 238.10	MS 238.10	10 Hole (L/R)
MT 238.11	MS 238.11	11 Hole (L/R)
MT 238.12	MS 238.12	12 Hole (L/R)

Meditech




Distal Medial Tibia Locking Plate

Primary Use : Lower and Tibia Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 239.03	MS 239.03	03 Hole (L/R)
MT 239.04	MS 239.04	04 Hole (L/R)
MT 239.05	MS 239.05	05 Hole (L/R)
MT 239.06	MS 239.06	06 Hole (L/R)
MT 239.07	MS 239.07	07 Hole (L/R)
MT 239.08	MS 239.08	08 Hole (L/R)
MT 239.09	MS 239.09	09 Hole (L/R)
MT 239.10	MS 239.10	10 Hole (L/R)
MT 239.11	MS 239.11	11 Hole (L/R)
MT 239.12	MS 239.12	12 Hole (L/R)
MT 239.13	MS 239.13	13 Hole (L/R)



Meditech



Low Band Distal Tibia Locking Plate


Primary Use : Lower and Tibia Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 240.03	MS 240.03	03 Hole (L/R)
MT 240.04	MS 240.04	04 Hole (L/R)
MT 240.05	MS 240.05	05 Hole (L/R)
MT 240.06	MS 240.06	06 Hole (L/R)
MT 240.07	MS 240.07	07 Hole (L/R)
MT 240.08	MS 240.08	08 Hole (L/R)
MT 240.09	MS 240.09	09 Hole (L/R)
MT 240.10	MS 240.10	10 Hole (L/R)
MT 240.11	MS 240.11	11 Hole (L/R)
MT 240.12	MS 240.12	12 Hole (L/R)
MT 240.13	MS 240.13	13 Hole (L/R)

Meditech



Distal Femur Locking Plate


Primary Use : Lower and Femur Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 241.03	MS 241.03	03 Hole (L/R)
MT 241.04	MS 241.04	04 Hole (L/R)
MT 241.05	MS 241.05	05 Hole (L/R)
MT 241.06	MS 241.06	06 Hole (L/R)
MT 241.07	MS 241.07	07 Hole (L/R)
MT 241.08	MS 241.08	08 Hole (L/R)
MT 241.09	MS 241.09	09 Hole (L/R)
MT 241.10	MS 241.10	10 Hole (L/R)
MT 241.11	MS 241.11	11 Hole (L/R)
MT 241.12	MS 241.12	12 Hole (L/R)
MT 241.13	MS 241.13	13 Hole (L/R)

Meditech



Antrolateral Distal Tibia Locking Plate


Primary Use : Lower and Tibia Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 242.03	MS 242.03	03 Hole (L/R)
MT 242.04	MS 242.04	04 Hole (L/R)
MT 242.05	MS 242.05	05 Hole (L/R)
MT 242.06	MS 242.06	06 Hole (L/R)
MT 242.07	MS 242.07	07 Hole (L/R)
MT 242.08	MS 242.08	08 Hole (L/R)
MT 242.09	MS 242.09	09 Hole (L/R)
MT 242.10	MS 242.10	10 Hole (L/R)
MT 242.11	MS 242.11	11 Hole (L/R)
MT 242.12	MS 242.12	12 Hole (L/R)
MT 242.13	MS 242.13	13 Hole (L/R)

Meditech



Philos Locking Plate Short


Primary Use : Proximal Humerus Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 243.03	MS 243.03	03 Hole
MT 243.04	MS 243.04	04 Hole
MT 243.05	MS 243.05	05 Hole
MT 243.06	MS 243.06	06 Hole
MT 243.07	MS 243.07	07 Hole
MT 243.08	MS 243.08	08 Hole
MT 243.09	MS 243.09	09 Hole
MT 243.10	MS 243.10	10 Hole
MT 243.11	MS 243.11	11 Hole
MT 243.12	MS 243.12	12 Hole
MT 243.13	MS 243.13	13 Hole

Meditech



Philos Locking Plate Long


Primary Use : Proximal Humerus Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 244.03	MS 244.03	03 Hole
MT 244.04	MS 244.04	04 Hole
MT 244.05	MS 244.05	05 Hole
MT 244.06	MS 244.06	06 Hole
MT 244.07	MS 244.07	07 Hole
MT 244.08	MS 244.08	08 Hole
MT 244.09	MS 244.09	09 Hole
MT 244.10	MS 244.10	10 Hole
MT 244.11	MS 244.11	11 Hole
MT 244.12	MS 244.12	12 Hole
MT 244.13	MS 244.13	13 Hole

Meditech



Proximal Humerus Locking Plate


Primary Use : Proximal Humerus Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 245.03	MS 245.03	03 Hole
MT 245.04	MS 245.04	04 Hole
MT 245.05	MS 245.05	05 Hole
MT 245.06	MS 245.06	06 Hole
MT 245.07	MS 245.07	07 Hole
MT 245.08	MS 245.08	08 Hole
MT 245.09	MS 245.09	09 Hole
MT 245.10	MS 245.10	10 Hole
MT 245.11	MS 245.11	11 Hole
MT 245.12	MS 245.12	12 Hole
MT 245.13	MS 245.13	13 Hole

Meditech



Cloverleaf Locking Plate


Primary Use : Proximal Humerus and Distal Tibia Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 246.03	MS 246.03	03 Hole (L/R)
MT 246.04	MS 246.04	04 Hole (L/R)
MT 246.05	MS 246.05	05 Hole (L/R)
MT 246.06	MS 246.06	06 Hole (L/R)
MT 246.07	MS 246.07	07 Hole (L/R)
MT 246.08	MS 246.08	08 Hole (L/R)
MT 246.09	MS 246.09	09 Hole (L/R)
MT 246.10	MS 246.10	10 Hole (L/R)
MT 246.11	MS 246.11	11 Hole (L/R)
MT 246.12	MS 246.12	12 Hole (L/R)
MT 246.13	MS 246.13	13 Hole (L/R)

Meditech



Clavicle Medial Locking Plate


Primary Use : Clavicle Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw
4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 247.05	MS 247.05	05 Hole (L/R)
MT 247.06	MS 247.06	06 Hole (L/R)
MT 247.07	MS 247.07	07 Hole (L/R)
MT 247.08	MS 247.08	08 Hole (L/R)
MT 247.09	MS 247.09	09 Hole (L/R)
MT 247.10	MS 247.10	10 Hole (L/R)



Meditech



Clavicle Hook Locking Plate

Primary Use : Clavicle Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw
4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 248.02	MS 248.02	02 Hole (L/R)
MT 248.03	MS 248.03	03 Hole (L/R)
MT 248.04	MS 248.04	04 Hole (L/R)
MT 248.05	MS 248.05	05 Hole (L/R)
MT 248.06	MS 248.06	06 Hole (L/R)
MT 248.07	MS 248.07	07 Hole (L/R)



Meditech



Superior Clavicle Locking Plate


Primary Use : Clavicle Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw 2.7mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 249.02	MS 249.02	02 Hole (L/R)
MT 249.03	MS 249.03	03 Hole (L/R)
MT 249.04	MS 249.04	04 Hole (L/R)
MT 249.05	MS 249.05	05 Hole (L/R)
MT 249.06	MS 249.06	06 Hole (L/R)
MT 249.07	MS 249.07	07 Hole (L/R)



Meditech



Y-Adaption Locking Plate

Primary Use : Distal Humerus Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM
(Prod. Code)

SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 250.02	MS 250.02	02 Hole (L/R)
MT 250.03	MS 250.03	03 Hole (L/R)
MT 250.04	MS 250.04	04 Hole (L/R)
MT 250.05	MS 250.05	05 Hole (L/R)
MT 250.06	MS 250.06	06 Hole (L/R)
MT 250.07	MS 250.07	07 Hole (L/R)



Meditech



Calcaneal Locking Plate

Primary Use : Calcaneal Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM
(Prod. Code)

SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 251.60

MS 251.60

60mm (L/R)


MT 251.70

MS 251.70

70mm (L/R)



Meditech




3.5mm Pilon Locking Plate

Primary Use : _____

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw
4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 252.02	MS 252.02	02 Hole (L/R)
MT 252.03	MS 252.03	03 Hole (L/R)
MT 252.04	MS 252.04	04 Hole (L/R)
MT 252.05	MS 252.05	05 Hole (L/R)
MT 252.06	MS 252.06	06 Hole (L/R)
MT 252.07	MS 252.07	07 Hole (L/R)
MT 252.08	MS 252.08	08 Hole (L/R)
MT 252.09	MS 252.09	09 Hole (L/R)
MT 252.10	MS 252.10	10 Hole (L/R)
MT 252.11	MS 252.11	11 Hole (L/R)
MT 252.12	MS 252.12	12 Hole (L/R)
MT 252.13	MS 252.13	13 Hole (L/R)

Meditech




3.5mm Pediatric Locking Plate

Primary Use : Pediatric Proximal Femur Fractures

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 253.03	MS 253.03	03 Hole (L/R)
MT 253.04	MS 253.04	04 Hole (L/R)
MT 253.05	MS 253.05	05 Hole (L/R)
MT 253.06	MS 253.06	06 Hole (L/R)
MT 253.07	MS 253.07	07 Hole (L/R)
MT 253.08	MS 253.08	08 Hole (L/R)
MT 253.09	MS 253.09	09 Hole (L/R)
MT 253.10	MS 253.10	10 Hole (L/R)
MT 253.11	MS 253.11	11 Hole (L/R)
MT 253.12	MS 253.12	12 Hole (L/R)
MT 253.13	MS 253.13	13 Hole (L/R)



Meditech



Pediatric Hip Locking Plate

Primary Use : Pediatric Hip Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 254.03	MS 254.03	03 Hole (L/R)
MT 254.04	MS 254.04	04 Hole (L/R)
MT 254.05	MS 254.05	05 Hole (L/R)
MT 254.06	MS 254.06	06 Hole (L/R)
MT 254.07	MS 254.07	07 Hole (L/R)
MT 254.08	MS 254.08	08 Hole (L/R)
MT 254.09	MS 254.09	09 Hole (L/R)
MT 254.10	MS 254.10	10 Hole (L/R)
MT 254.11	MS 254.11	11 Hole (L/R)
MT 254.12	MS 254.12	12 Hole (L/R)
MT 254.13	MS 254.13	13 Hole (L/R)

Meditech



Posterior Medial Tibia Locking Plate


Primary Use : Posterior Medial Tibia Fracture

Matching Screw : 3.5mm Cortical and 3.5mm Locking Screw

4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 255.03	MS 255.03	03 Hole (L/R)
MT 255.04	MS 255.04	04 Hole (L/R)
MT 255.05	MS 255.05	05 Hole (L/R)
MT 255.06	MS 255.06	06 Hole (L/R)
MT 255.07	MS 255.07	07 Hole (L/R)
MT 255.08	MS 255.08	08 Hole (L/R)
MT 255.09	MS 255.09	09 Hole (L/R)
MT 255.10	MS 255.10	10 Hole (L/R)
MT 255.11	MS 255.11	11 Hole (L/R)
MT 255.12	MS 255.12	12 Hole (L/R)
MT 255.13	MS 255.13	13 Hole (L/R)

Meditech




4.5mm Recon Locking Plate


Primary Use : Pelvic Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 256.04	MS 256.04	04 Hole
MT 256.05	MS 256.05	05 Hole
MT 256.06	MS 256.06	06 Hole
MT 256.07	MS 256.07	07 Hole
MT 256.08	MS 256.08	08 Hole
MT 256.09	MS 256.09	09 Hole
MT 256.10	MS 256.10	10 Hole
MT 256.11	MS 256.11	11 Hole
MT 256.12	MS 256.12	12 Hole
MT 256.13	MS 256.13	13 Hole



Meditech




4.5mm Narrow Locking Plate


Primary Use : Humerus Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 257.04	MS 257.04	04 Hole
MT 257.05	MS 257.05	05 Hole
MT 257.06	MS 257.06	06 Hole
MT 257.07	MS 257.07	07 Hole
MT 257.08	MS 257.08	08 Hole
MT 257.09	MS 257.09	09 Hole
MT 257.10	MS 257.10	10 Hole
MT 257.11	MS 257.11	11 Hole
MT 257.12	MS 257.12	12 Hole
MT 257.13	MS 257.13	13 Hole



Meditech




4.5mm Broad Locking Plate Straight

Primary Use : Femur Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 258.04	MS 258.04	04 Hole
MT 258.05	MS 258.05	05 Hole
MT 258.06	MS 258.06	06 Hole
MT 258.07	MS 258.07	07 Hole
MT 258.08	MS 258.08	08 Hole
MT 258.09	MS 258.09	09 Hole
MT 258.10	MS 258.10	10 Hole
MT 258.11	MS 258.11	11 Hole
MT 258.12	MS 258.12	12 Hole
MT 258.13	MS 258.13	13 Hole
MT 258.14	MS 258.14	14 Hole
MT 258.14	MS 258.15	15 Hole
MT 258.16	MS 258.16	16 Hole




4.5mm Broad Locking Plate Curved

Primary Use : Femur Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 259.04	MS 259.04	04 Hole
MT 259.05	MS 259.05	05 Hole
MT 259.06	MS 259.06	06 Hole
MT 259.07	MS 259.07	07 Hole
MT 259.08	MS 259.08	08 Hole
MT 259.09	MS 259.09	09 Hole
MT 259.10	MS 259.10	10 Hole
MT 259.11	MS 259.11	11 Hole
MT 259.12	MS 259.12	12 Hole
MT 259.13	MS 259.13	13 Hole
MT 259.14	MS 259.14	14 Hole
MT 259.15	MS 259.15	15 Hole
MT 259.16	MS 259.16	16 Hole
MT 259.17	MS 259.17	17 Hole
MT 259.18	MS 259.18	18 Hole




Lateral Tibia Locking Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 260.03	MS 260.03	03 Hole (L/R)
MT 260.04	MS 260.04	04 Hole (L/R)
MT 260.05	MS 260.05	05 Hole (L/R)
MT 260.06	MS 260.06	06 Hole (L/R)
MT 260.07	MS 260.07	07 Hole (L/R)
MT 260.08	MS 260.08	08 Hole (L/R)
MT 260.09	MS 260.09	09 Hole (L/R)
MT 260.10	MS 260.10	10 Hole (L/R)
MT 260.11	MS 260.11	11 Hole (L/R)
MT 260.12	MS 260.12	12 Hole (L/R)
MT 260.13	MS 260.13	13 Hole (L/R)

Meditech




Proximal Tibia Locking Raft Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 4.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 261.03	MS 261.03	03 Hole (L/R)
MT 261.04	MS 261.04	04 Hole (L/R)
MT 261.05	MS 261.05	05 Hole (L/R)
MT 261.06	MS 261.06	06 Hole (L/R)
MT 261.07	MS 261.07	07 Hole (L/R)
MT 261.08	MS 261.08	08 Hole (L/R)
MT 261.09	MS 261.09	09 Hole (L/R)
MT 261.10	MS 261.10	10 Hole (L/R)
MT 261.11	MS 261.11	11 Hole (L/R)
MT 261.12	MS 261.12	12 Hole (L/R)
MT 261.13	MS 261.13	13 Hole (L/R)

Meditech




T-Buttress Locking Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 262.03	MS 262.03	03 Hole
MT 262.04	MS 262.04	04 Hole
MT 262.05	MS 262.05	05 Hole
MT 262.06	MS 262.06	06 Hole
MT 262.07	MS 262.07	07 Hole
MT 262.08	MS 262.08	08 Hole
MT 262.09	MS 262.09	09 Hole
MT 262.10	MS 262.10	10 Hole
MT 262.11	MS 262.11	11 Hole
MT 262.12	MS 262.12	12 Hole
MT 262.13	MS 262.13	13 Hole

Meditech




Tomofix Medial Tibia Locking Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 263.03	MS 263.03	03 Hole
MT 263.04	MS 263.04	04 Hole
MT 263.05	MS 263.05	05 Hole
MT 263.06	MS 263.06	06 Hole
MT 263.07	MS 263.07	07 Hole
MT 263.08	MS 263.08	08 Hole
MT 263.09	MS 263.09	09 Hole
MT 263.10	MS 263.10	10 Hole
MT 263.11	MS 263.11	11 Hole
MT 263.12	MS 263.12	12 Hole
MT 263.13	MS 263.13	13 Hole

Meditech



Osteotomi Medial Distal Femur Locking Plate

Primary Use : Distal Femur Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM
(Prod. Code)

SS. 316L
Prod. Code

SPECIFICATION
(Holes)

MT 264.03	MS 264.03	03 Hole (L/R)
MT 264.04	MS 264.04	04 Hole (L/R)
MT 264.05	MS 264.05	05 Hole (L/R)
MT 264.06	MS 264.06	06 Hole (L/R)
MT 264.07	MS 264.07	07 Hole (L/R)
MT 264.08	MS 264.08	08 Hole (L/R)
MT 264.09	MS 264.09	09 Hole (L/R)
MT 264.10	MS 264.10	10 Hole (L/R)
MT 264.11	MS 264.11	11 Hole (L/R)
MT 264.12	MS 264.12	12 Hole (L/R)
MT 264.13	MS 264.13	13 Hole (L/R)

Meditech




L-Buttress Locking Plate**Primary Use :** Upper Tibia Fracture**Matching Screw :** 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

65

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 265.03	MS 265.03	03 Hole (L/R)
MT 265.04	MS 265.04	04 Hole (L/R)
MT 265.05	MS 265.05	05 Hole (L/R)
MT 265.06	MS 265.06	06 Hole (L/R)
MT 265.07	MS 265.07	07 Hole (L/R)
MT 265.08	MS 265.08	08 Hole (L/R)
MT 265.09	MS 265.09	09 Hole (L/R)
MT 265.10	MS 265.10	10 Hole (L/R)
MT 265.11	MS 265.11	11 Hole (L/R)
MT 265.12	MS 265.12	12 Hole (L/R)
MT 265.13	MS 265.13	13 Hole (L/R)

Meditech



L-Locking Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 266.03	MS 266.03	03 Hole (L/R)
MT 266.04	MS 266.04	04 Hole (L/R)
MT 266.05	MS 266.05	05 Hole (L/R)
MT 266.06	MS 266.06	06 Hole (L/R)
MT 266.07	MS 266.07	07 Hole (L/R)
MT 266.08	MS 266.08	08 Hole (L/R)
MT 266.09	MS 266.09	09 Hole (L/R)
MT 266.10	MS 266.10	10 Hole (L/R)
MT 266.11	MS 266.11	11 Hole (L/R)
MT 266.12	MS 266.12	12 Hole (L/R)
MT 266.13	MS 266.13	13 Hole (L/R)




T-Locking Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 267.03	MS 267.03	03 Hole
MT 267.04	MS 267.04	04 Hole
MT 267.05	MS 267.05	05 Hole
MT 267.06	MS 267.06	06 Hole
MT 267.07	MS 267.07	07 Hole
MT 267.08	MS 267.08	08 Hole
MT 267.09	MS 267.09	09 Hole
MT 267.10	MS 267.10	10 Hole
MT 267.11	MS 267.11	11 Hole
MT 267.12	MS 267.12	12 Hole
MT 267.13	MS 267.13	13 Hole

Meditech



Proximal Tibia L-Locking Plate

Primary Use : Upper Tibia Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 268.03	MS 268.03	03 Hole (L/R)
MT 268.04	MS 268.04	04 Hole (L/R)
MT 268.05	MS 268.05	05 Hole (L/R)
MT 268.06	MS 268.06	06 Hole (L/R)
MT 268.07	MS 268.07	07 Hole (L/R)
MT 268.08	MS 268.08	08 Hole (L/R)
MT 268.09	MS 268.09	09 Hole (L/R)
MT 268.10	MS 268.10	10 Hole (L/R)
MT 268.11	MS 268.11	11 Hole (L/R)
MT 268.12	MS 268.12	12 Hole (L/R)
MT 268.13	MS 268.13	13 Hole (L/R)


Meditech



Distal Femur Locking Plate**Primary Use :** Distal Femur Fracture**Matching Screw :** 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

69

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 269.03	MS 269.03	03 Hole (L/R)
MT 269.04	MS 269.04	04 Hole (L/R)
MT 269.05	MS 269.05	05 Hole (L/R)
MT 269.06	MS 269.06	06 Hole (L/R)
MT 269.07	MS 269.07	07 Hole (L/R)
MT 269.08	MS 269.08	08 Hole (L/R)
MT 269.09	MS 269.09	09 Hole (L/R)
MT 269.10	MS 269.10	10 Hole (L/R)
MT 269.11	MS 269.11	11 Hole (L/R)
MT 269.12	MS 269.12	12 Hole (L/R)
MT 269.13	MS 269.13	13 Hole (L/R)
MT 269.14	MS 269.14	14 Hole (L/R)
MT 269.15	MS 269.15	15 Hole (L/R)




Distal Femur Medial Locking Plate

Primary Use : Distal Femur Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 270.03	MS 270.03	03 Hole (L/R)
MT 270.04	MS 270.04	04 Hole (L/R)
MT 270.05	MS 270.05	05 Hole (L/R)
MT 270.06	MS 270.06	06 Hole (L/R)
MT 270.07	MS 270.07	07 Hole (L/R)
MT 270.08	MS 270.08	08 Hole (L/R)
MT 270.09	MS 270.09	09 Hole (L/R)
MT 270.10	MS 270.10	10 Hole (L/R)
MT 270.11	MS 270.11	11 Hole (L/R)
MT 270.12	MS 270.12	12 Hole (L/R)
MT 270.13	MS 270.13	13 Hole (L/R)


Meditech



Condylar Buttress Locking Plate**Primary Use :** Condylar Femur Fracture**Matching Screw :** 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

71

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 271.03	MS 271.03	03 Hole (L/R)
MT 271.04	MS 271.04	04 Hole (L/R)
MT 271.05	MS 271.05	05 Hole (L/R)
MT 271.06	MS 271.06	06 Hole (L/R)
MT 271.07	MS 271.07	07 Hole (L/R)
MT 271.08	MS 271.08	08 Hole (L/R)
MT 271.09	MS 271.09	09 Hole (L/R)
MT 271.10	MS 271.10	10 Hole (L/R)
MT 271.11	MS 271.11	11 Hole (L/R)
MT 271.12	MS 271.12	12 Hole (L/R)
MT 271.13	MS 271.13	13 Hole (L/R)
MT 271.14	MS 271.14	14 Hole (L/R)
MT 271.15	MS 271.15	15 Hole (L/R)




95°DCS Locking Plate

Primary Use : Condylar Femur Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw and 12mm leg Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 272.03	MS 272.03	03 Hole
MT 272.04	MS 272.04	04 Hole
MT 272.05	MS 272.05	05 Hole
MT 272.06	MS 272.06	06 Hole
MT 272.07	MS 272.07	07 Hole
MT 272.08	MS 272.08	08 Hole
MT 272.09	MS 272.09	09 Hole
MT 272.10	MS 272.10	10 Hole
MT 272.11	MS 272.11	11 Hole
MT 272.12	MS 272.12	12 Hole
MT 272.13	MS 272.13	13 Hole
MT 272.14	MS 272.14	14 Hole
MT 272.15	MS 272.15	15 Hole




DHS Locking Plate 125°, 130, 135 & 140 °

Primary Use : Hip Fractures

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw and 12mm leg Screw

73

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 273.03	MS 273.03	03 Hole
MT 273.04	MS 273.04	04 Hole
MT 273.05	MS 273.05	05 Hole
MT 273.06	MS 273.06	06 Hole
MT 273.07	MS 273.07	07 Hole
MT 273.08	MS 273.08	08 Hole
MT 273.09	MS 273.09	09 Hole
MT 273.10	MS 273.10	10 Hole
MT 273.11	MS 273.11	11 Hole
MT 273.12	MS 273.12	12 Hole
MT 273.13	MS 273.13	13 Hole
MT 273.14	MS 273.14	14 Hole
MT 273.15	MS 273.15	15 Hole



Trochanteric Stabilizing Plate**Primary Use :** Trochanteric Hip Fractures**Matching Screw :** 4.5mm Cortical and 5.0mm Locking Screw and 12mm leg Screw

74

TITANIUM
(Prod. Code)**SS. 316L**
Prod. Code**SPECIFICATION**
(Holes)


MT 274.03

MS 274.03

03 Hole



Meditech




Proximal Femur Hook Locking Plate

Primary Use : Proximal Femur Neck Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 6.5mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 275.03	MS 275.03	03 Hole (L/R)
MT 275.04	MS 275.04	04 Hole (L/R)
MT 275.05	MS 275.05	05 Hole (L/R)
MT 275.06	MS 275.06	06 Hole (L/R)
MT 275.07	MS 275.07	07 Hole (L/R)
MT 275.08	MS 275.08	08 Hole (L/R)
MT 275.09	MS 275.09	09 Hole (L/R)
MT 275.10	MS 275.10	10 Hole (L/R)
MT 275.11	MS 275.11	11 Hole (L/R)
MT 275.12	MS 275.12	12 Hole (L/R)
MT 275.13	MS 275.13	13 Hole (L/R)

Meditech




Proximal Femur Lateral Locking Plate

Primary Use : Proximal Femur Lateral Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 5.0mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 276.03	MS 276.03	03 Hole (L/R)
MT 276.04	MS 276.04	04 Hole (L/R)
MT 276.05	MS 276.05	05 Hole (L/R)
MT 276.06	MS 276.06	06 Hole (L/R)
MT 276.07	MS 276.07	07 Hole (L/R)
MT 276.08	MS 276.08	08 Hole (L/R)
MT 276.09	MS 276.09	09 Hole (L/R)
MT 276.10	MS 276.10	10 Hole (L/R)
MT 276.11	MS 276.11	11 Hole (L/R)
MT 276.12	MS 276.12	12 Hole (L/R)
MT 276.13	MS 276.13	13 Hole (L/R)

Meditech




Proximal Femur Locking Plate

Primary Use : Proximal Femur Fracture

Matching Screw : 4.5mm Cortical and 5.0mm Locking Screw
and 6.5mm Cancellous Locking Screw

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 277.03	MS 277.03	03 Hole (L/R)
MT 277.04	MS 277.04	04 Hole (L/R)
MT 277.05	MS 277.05	05 Hole (L/R)
MT 277.06	MS 277.06	06 Hole (L/R)
MT 277.07	MS 277.07	07 Hole (L/R)
MT 277.08	MS 277.08	08 Hole (L/R)
MT 277.09	MS 277.09	09 Hole (L/R)
MT 277.10	MS 277.10	10 Hole (L/R)
MT 277.11	MS 277.11	11 Hole (L/R)
MT 277.12	MS 277.12	12 Hole (L/R)
MT 277.13	MS 277.13	13 Hole (L/R)

Meditech



2.0mm Locking Head Screw

Primary Use : Mini Fracture Plate

Drill Bit Used : 1.5mm Drill Bit

78

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 278.06	MS 278.06	06mm
MT 278.08	MS 278.08	08mm
MT 278.10	MS 278.10	10mm
MT 278.12	MS 278.12	12mm
MT 278.14	MS 278.14	14mm
MT 278.16	MS 278.16	16mm
MT 278.18	MS 278.18	18mm
MT 278.20	MS 278.20	20mm
MT 278.22	MS 278.22	22mm
MT 278.24	MS 278.24	24mm
MT 278.26	MS 278.26	26mm
MT 278.28	MS 278.28	28mm
MT 278.30	MS 278.30	30mm
MT 278.35	MS 278.35	35mm
MT 278.40	MS 278.40	40mm



2.7mm Locking Head Screw**Primary Use :** Mini Fracture Plate**Drill Bit Used :** 2.0mm Drill Bit

79

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
---------------------------------	-------------------------------	---------------------------------

MT 279.08	MS 279.08	08mm
MT 279.10	MS 279.10	10mm
MT 279.12	MS 279.12	12mm
MT 279.14	MS 279.14	14mm
MT 279.16	S 279.16	16mm
MT 279.18	MS 279.18	18mm
MT 279.20	MS 279.20	20mm
MT 279.22	MS 279.22	22mm
MT 279.24	MS 279.24	24mm
MT 279.26	MS 279.26	26mm
MT 279.28	MS 279.28	28mm
MT 279.30	MS 279.30	30mm
MT 279.32	MS 279.32	32mm
MT 279.34	MS 279.34	34mm
MT 279.36	MS 279.36	36mm
MT 279.38	MS 279.38	38mm
MT 279.40	MS 279.40	40mm
MT 279.45	MS 279.45	45mm
MT 279.50	MS 279.50	50mm



3.5mm Locking Head Screw

Primary Use : Small Fracture Plate

Drill Bit Used : 2.7mm Drill Bit

80

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 280.12	MS 280.12	12mm
MT 280.14	MS 280.14	14mm
MT 280.16	MS 280.16	16mm
MT 280.18	MS 280.18	18mm
MT 280.20	MS 280.20	20mm
MT 280.22	MS 280.22	22mm
MT 280.24	MS 280.24	24mm
MT 280.26	MS 280.26	26mm
MT 280.28	MS 280.28	28mm
MT 280.30	MS 280.30	30mm
MT 280.32	MS 280.32	32mm
MT 280.34	MS 280.34	34mm
MT 280.36	MS 280.36	36mm
MT 280.38	MS 280.38	38mm
MT 280.40	MS 280.40	40mm
MT 280.45	MS 280.45	45mm
MT 280.50	MS 280.50	50mm
MT 280.55	MS 280.55	55mm
MT 280.60	MS 280.60	60mm
MT 280.65	MS 280.65	65mm
MT 280.70	MS 280.70	70mm



5.0mm Locking Head Screw


Primary Use : Large Fracture Plate

Drill Bit Used : 4.0mm Drill Bit

81

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
--------------------------	------------------------	--------------------------

MT 281.18	MS 281.18	18mm
MT 281.20	MS 281.20	20mm
MT 281.22	MS 281.22	22mm
MT 281.24	MS 281.24	24mm
MT 281.26	MS 281.26	26mm
MT 281.28	MS 281.28	28mm
MT 281.30	MS 281.30	30mm
MT 281.32	MS 281.32	32mm
MT 281.34	MS 281.34	34mm
MT 281.36	MS 281.36	36mm
MT 281.38	MS 281.38	38mm
MT 281.40	MS 281.40	40mm
MT 281.42	MS 281.42	42mm
MT 281.44	MS 281.44	44mm
MT 281.46	MS 281.46	46mm
MT 281.48	MS 281.48	48mm
MT 281.50	MS 281.50	50mm
MT 281.55	MS 281.55	55mm
MT 281.60	MS 281.60	60mm
MT 281.65	MS 281.65	65mm
MT 281.70	MS 281.70	70mm
MT 281.75	MS 281.75	75mm
MT 281.80	MS 281.80	80mm



2.7mm Cancellous Locking Head Screw

Primary Use : Mini Fracture Plate

Drill Bit Used : 2.0mm Drill Bit

82

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
--------------------------	------------------------	--------------------------

MT 282.08	MS 282.08	08mm
MT 282.10	MS 282.10	10mm
MT 282.12	MS 282.12	12mm
MT 282.14	MS 282.14	14mm
MT 282.16	MS 282.16	16mm
MT 282.18	MS 282.18	18mm
MT 282.20	MS 282.20	20mm
MT 282.22	MS 282.22	22mm
MT 282.24	MS 282.24	24mm
MT 282.26	MS 282.26	26mm
MT 282.28	MS 282.28	28mm
MT 282.30	MS 282.30	30mm
MT 282.32	MS 282.32	32mm
MT 282.34	MS 282.34	34mm
MT 282.36	MS 282.36	36mm
MT 282.38	MS 282.38	38mm
MT 282.40	MS 282.40	40mm
MT 282.45	MS 282.45	45mm
MT 282.50	MS 282.50	50mm




4.0mm Cancellous Locking Head Screw**Primary Use :** Small Fracture Plate**Drill Bit Used :** 2.7mm Drill Bit

83

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
---------------------------------	-------------------------------	---------------------------------

MT 283.12	MS 283.12	12mm
MT 283.14	MS 283.14	14mm
MT 283.16	MS 283.16	16mm
MT 283.18	MS 283.18	18mm
MT 283.20	MS 283.20	20mm
MT 283.22	MS 283.22	22mm
MT 283.24	MS 283.24	24mm
MT 283.26	MS 283.26	26mm
MT 283.28	MS 283.28	28mm
MT 283.30	MS 283.30	30mm
MT 283.35	MS 283.35	35mm
MT 283.40	MS 283.40	40mm
MT 283.45	MS 283.45	45mm
MT 283.50	MS 283.50	50mm
MT 283.55	MS 283.55	55mm
MT 283.60	MS 283.60	60mm
MT 283.65	MS 283.65	65mm
MT 283.70	MS 283.70	70mm
MT 283.75	MS 283.75	75mm
MT 283.80	MS 283.80	80mm



5.0mm Cancellous Locking Head Screw**Primary Use :** Large Fracture Plate**Drill Bit Used :** 4.0mm Drill Bit

84

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
MT 284.30	MS 284.30	30mm
MT 284.35	MS 284.35	35mm
MT 284.40	MS 284.40	40mm
MT 284.45	MS 284.45	45mm
MT 284.50	MS 284.50	50mm
MT 284.55	MS 284.55	55mm
MT 284.60	MS 284.60	60mm
MT 284.65	MS 284.65	65mm
MT 284.70	MS 284.70	70mm
MT 284.75	MS 284.75	75mm
MT 284.80	MS 284.80	80mm
MT 284.85	MS 284.85	85mm
MT 284.90	MS 284.90	90mm
MT 284.95	MS 284.95	95mm
MT 284.100	MS 284.100	100mm
MT 284.105	MS 284.105	105mm
MT 284.110	MS 284.110	110mm



6.5mm Cancellous Locking Head Screw**Primary Use :** Large Fracture Plate**Drill Bit Used :** 4.3mm Drill Bit

85

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
---------------------------------	-------------------------------	---------------------------------

MT 285.30	MS 285.30	30mm
MT 285.35	MS 285.35	35mm
MT 285.40	MS 285.40	40mm
MT 285.45	MS 285.45	45mm
MT 285.50	MS 285.50	50mm
MT 285.55	MS 285.55	55mm
MT 285.60	MS 285.60	60mm
MT 285.65	MS 285.65	65mm
MT 285.70	MS 285.70	70mm
MT 285.75	MS 285.75	75mm
MT 285.80	MS 285.80	80mm
MT 285.85	MS 285.85	85mm
MT 285.90	MS 285.90	90mm
MT 285.95	MS 285.95	95mm
MT 285.100	MS 285.100	100mm
MT 285.105	MS 285.105	105mm
MT 285.110	MS 285.110	110mm




7.3mm Cancellous Cannulated Locking Head Screw**Primary Use :** Large Fracture Plate**Drill Bit Used :** 4.3mm Cannulated Drill Bit

86

TITANIUM (Prod. Code)	SS. 316L Prod. Code	SPECIFICATION (Holes)
---------------------------------	-------------------------------	---------------------------------

MT 286.30	MS 286.30	30mm
MT 286.35	MS 286.35	35mm
MT 286.40	MS 286.40	40mm
MT 286.45	MS 286.45	45mm
MT 286.50	MS 286.50	50mm
MT 286.55	MS 286.55	55mm
MT 286.60	MS 286.60	60mm
MT 286.65	MS 286.65	65mm
MT 286.70	MS 286.70	70mm
MT 286.75	MS 286.75	75mm
MT 286.80	MS 286.80	80mm
MT 286.85	MS 286.85	85mm
MT 286.90	MS 286.90	90mm
MT 286.95	MS 286.95	95mm
MT 286.100	MS 286.100	100mm
MT 286.105	MS 286.105	105mm
MT 286.110	MS 286.110	110mm





Meditech

Meditech Ortho India Company



Meditech Ortho India Company

📍 Khasra No 28/19, House No 742, Ground Floor, Gali No - 5/12,
Village Jharoda Majra, Sangam Vihar, North East Delhi - 110084

☎ + 91-9716713263